# Gebruiksaanweisung CARECATH<sup>live</sup>, CARECATH<sup>leiv</sup>

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## GB Instruction for use

Hydrophilic coated, ready to use, sterile single use catheter of medical grade PVC for Intermittent catheterization.

### Specifications

The product is offered in the size of 06Fr, 08Fr, 10Fr, 12Fr, 14Fr, 16Fr and 18Fr, based on the Outer Diameter, and in the following lengths:

| Ref No.    | CH-Size          | Length |
|------------|------------------|--------|
| LIV-H0817  | 08Fr/CH   2,67mm | 17cm   |
| LIV-H1017  | 10Fr/CH   3,33mm | 17cm   |
| LIV-H1217  | 12Fr/CH   4mm    | 17cm   |
| LIV-H1417  | 14Fr/CH   4,67mm | 17cm   |
| LIV-H1617  | 16Fr/CH   5,33mm | 17cm   |
| LIV-H1817  | 18Fr/CH   6mm    | 17cm   |
| LEIV-H0641 | 06Fr/CH   2mm    | 41cm   |
| LEIV-H0841 | 08Fr/CH   2,67mm | 41cm   |
| LEIV-H1041 | 10Fr/CH   3,33mm | 41cm   |
| LEIV-H1241 | 12Fr/CH   4mm    | 41cm   |
| LEIV-H1441 | 14Fr/CH   4,67mm | 41cm   |
| LEIV-H1641 | 16Fr/CH   5,34mm | 41cm   |
| LEIV-H1841 | 18Fr   CH   6mm  | 41cm   |

All specifications are fully complied with ISO 20696:2018 standards.

### Indications for Use

Intermittent urinary catheterization by inserting through the urethra to pass urine from the bladder.

### Potential adverse events

- Bleeding in the urethra
- Catheter-associated urinary tract infection
- Urethral injuries and inflammation

### Contraindications

- Acute urethritis
- Acute prostatitis
- Acute epididymitis
- Patients with PVC allergy
- Patients have calcareous urolithiasis

### Precautions

- This product can be used in healthcare centers under supervision of appropriately qualified and trained personnel and/or at home by patients according to the instructions for use.
- This product is disposable, do not reuse.
- To prevent urinary tract infection, appropriate aseptic steps should be taken to areas around the urethral opening prior to insertion.
- Do not use if the package is damaged, contaminated, damaged or expired.
- Special attention shall be taken for patients with disturbance of consciousness as involuntary pull

- of the catheter may lead to bladder and urethra mucosal injury or break of catheter.
- Intubation and extubation should be as gentle as possible to prevent urethral injury.
- Slowly advance the catheter through the urethra into the bladder. If substantial resistance is met, do not force the catheter. Slowly pull out the catheter and look for reasons.
- If there is no urine draining into your bag, or toilet, take the following steps: Check for and remove any kinks in the catheter or the drainage bag tubing. Check the position of your catheter and drainage bag. Ensure the bag is positioned below your bladder when you are lying, sitting or standing.
- For overfilled bladder, the urine should not be discharged too fast which may lead to shock or bladder bleeding.

### Warning

After infiltrating the catheter, do not wipe or scrub the surface of the catheter by either physical force or chemicals to avoid damaging the hydrophilic coating which may result in efficacy and gliding ability loss.

This is a single use device. Do not re-sterilize any portion of this device. Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of the device, which may lead to device failure, and/or lead to injury, illness or death of the patient.

### Before use:

Select an appropriate specification of catheter which is suitable for the specific procedure. Carefully inspect the package and catheter prior to use. Discard and replace the product if it is damaged. The operator shall follow the standard aseptic procedure being applied to its own and its patient before conducting the procedure.

### Method of use:

- Wash your hands thoroughly and disinfect them.
  - Open the package and expose the funnel.
  - Holding the funnel, carefully pull the catheter out of the packaging.
  - Carefully insert the catheter through the urethral orifice and advance into the bladder until urine starts flow out after the eyelets of the catheter reach the bladder cavity.
  - If the patient is located in front of the toilet facility, the urine can drain into the toilet. If the patient is lying down, place the funnel end of the catheter into a container or connect a drainage bag by inserting the plug of drainage bag into the catheter funnel to collect the urine. Hold the catheter until urine stops flowing.
  - When urine stops flowing, begin to withdraw the catheter. It is recommended to slowly rotate the catheter during withdrawal, stopping each time urine begins to flow until the bladder is empty.
  - After use, remove the catheter completely from the urethra and dispose it in accordance with local administrative regulations to prevent environmental pollution and potential health hazards.
  - Wash hands after catheterization.
- \* 5, 6 steps are universal for both male and female patients.

### Validity period and storage method

- The product is valid for 3 years. For single use only, do not repeat sterilization and reuse.

- Use immediately after the package is opened and dispose after use.

- The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures.

### Sterilization method

- Product is sterilized by irradiation.

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