

Instructions for Use

**Name**

Disposable bib

Model

38*50cm, 40*60cm, 48*58cm, 36*64+10cm, 36*66+10cm, 37*50+10cm, 37*64+10cm, 37*65+10cm, 37*66+10cm, 38*65+10cm 38*70+10cm, 40*60+10cm, 45*60+10cm

Intended Use

The product is intended to protect the body and clothing from oral liquid, food, water, etc., polluting. This is a single use, disposable device(s), provided non-sterile.

Storage

The product should be stored in a cool dry area, away from heat and direct sunlight.

Use Instruction

- 1.Read the label and instruction for use carefully before used;
- 2.Take off the package, do not use the product once the package was damaged or the product has expired;
- 3.Place the disposable bib on the chest below the neck.

Shelf Life

5 years

Caution


1. It is forbidden to use when the packaging or the product is damaged;
2. This product is a single-use product and it is forbidden to reuse;
3. Replace the product in time if it is damaged or penetrated.

Disposal

Please dispose the product after use to comply with local regulation.

Labels, Packing Logo Design

Symbol	Introductions	Symbol	Introductions
	CE Symbol		medical device
	Manufacturer Name Address		AUTHORISED REPRESENTATIVE IN THE EUROPEAN COMMUNITY

	Manufacture Date		Use until year & month (Expiration date)
	Batch Code		Do not reuse" are "single use, "Use only once
	Warnings and Precautions		non-sterile
	Consult instructions for use		Do not use if package is damaged
	Country of Manufacture		Model Number
	Importer		Distributor
	Unique Device Identification		

Manufacturer Information



Company: Changzhou Newrise Medical & Hygiene Products Co., Ltd.
Address: No.39 Chuangye Road, Zhulin Town, Jintan District, Changzhou
City, Jiangsu, China 213241

European Authorized Representative



Company: SUNGO Europe B.V.
Address: Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan
den IJssel, The Netherlands