



**System provider  
for Incontinence Care**

**CATALOGUE**  
[www.ghc-medical.de](http://www.ghc-medical.de)

## Introduction

With our assortment of revulsive incontinence products, we as GHC German Health Care GmbH focus on an area, where our mission is to offer our customers an extensive range of different products as well as offering a very high level of expertise. And our customers can rely on that – Even during challenging times!

Within our product assortment, you will find a wide selection of transurethral and suprapubic balloon catheters, different catheter insertion kits, a wide range of urine collection systems and bags as well as bag holding systems and accessories.

We have positioned ourselves well for the future by having expanded our warehouse capacities, added a new warehouse management software as well as keeping up the high level of service orientation with various service offers. For example, we have created a customer portal where all relevant MDR-documents are available for you to view or download at any time. This service saves you a lot of time and gives you the benefit of always having access to all necessary documents.

In our current catalogue you will find new products that we have added to our assortment for our customers ranging from homecare companies, medical supply stores and clinics. For instance, we are introducing our new CARE *CON* series (page 38-40), which are Male External Catheters (MEC) in a very good and reliable quality at an attractive price. Further, we are introducing our new transurethral soft silicone balloon catheter (page 9), which due to the softer shaft material is an excellent alternative to the current balloon catheter range in the homecare market. Alongside that, we are implementing an additional new transurethral silicone catheter with a flute tip (page 14). This variant is a good solution for patients with problematic urine findings. Lastly, we also have implemented new suprapubic balloon catheters as alternatives (pages 25 and 29) to our previous product line-up. Besides the high product quality, we also wanted to enhance our delivery capabilities by adding the new products to our assortment.

With this catalogue we offer you an overview of the complete range of GHC's product portfolio. You are warmly invited to visit **[www.ghc-medical.de](http://www.ghc-medical.de)** for up-to-date information.

If our offer has sparked your interest, feel free to get in touch with us.

We look forward talking to you.

Warmest greetings from the GHC management,

  
Dirk Dohmann  
Kay Kruse

Table of Contents

|                                                                                              |          |
|----------------------------------------------------------------------------------------------|----------|
| Introduction                                                                                 | page 3   |
| Table of Contents                                                                            | page 4,5 |
| Quality policy of GHC                                                                        | page 6   |
| Catheter guide                                                                               | page 7   |
|                                                                                              |          |
| <b>CARE FLOW Transurethral catheter</b>                                                      |          |
| CARE FLOW 2-way 100% silicone catheter, nelaton tip                                          | page 8   |
| CARE FLOW SOFT- 2-way Soft 100% silicone catheter, nelaton tip <b>new</b>                    | page 9   |
| CARE FLOW 2-way integral 100% silicone catheter, nelaton tip and plug <b>new</b>             | page 10  |
| BLUECATH 2-way integral 100% silicone catheter, nelaton tip and plug <b>new</b>              | page 11  |
| CARE FLOW 2-way silicone elastomer catheter, nelaton tip                                     | page 12  |
| CARE FLOW 2-way 100% silicone catheter, tiemann tip <b>new</b>                               | page 13  |
| CARE FLOW 2-way 100% silicone catheter, flute tip <b>new</b>                                 | page 14  |
| CARE FLOW 2-way silicone elastomer catheter, couvelaire tip                                  | page 15  |
| CARE FLOW 2-way siliconised latex catheter, nelaton tip                                      | page 16  |
| CARE FLOW PVC single-use female catheter, nelaton tip                                        | page 17  |
| CARE FLOW PVC single-use male catheter, nelaton and tiemann tip                              | page 18  |
|                                                                                              |          |
| <b>CARE FLOW Irrigation catheter</b>                                                         |          |
| CARE FLOW 3-way 100% silicone irrigation catheter, nelaton tip                               | page 19  |
| CARE FLOW 3-way 100% silicone irrigation catheter, dufour tip                                | page 20  |
| CARE FLOW 3-way 100% silicone irrigation catheter, couvelaire tip                            | page 21  |
| CARE FLOW 3-way 100% silicone irrigation catheter, mercier tip                               | page 22  |
|                                                                                              |          |
| <b>CARE FLOW Suprapubic catheter</b>                                                         |          |
| CARE FLOW UNIBAL 2-way integral 100% silicone catheter, central open tip                     | page 23  |
| CARE FLOW UNIBAL 2-way integral 100% silicone catheter, central open tip, with plug          | page 24  |
| CARE FLOW 2-way 100% silicone catheter integral, central open tip, with plug <b>new</b>      | page 25  |
| CARE FLOW 2-way integral 100% silicone catheter, central open tip, short                     | page 26  |
| BLUECATH 2-way integral 100% silicone catheter, dilatation tip, with plug                    | page 27  |
| BLUECATH 2-way integral 100% silicone catheter, central open tip, with plug                  | page 28  |
| CARE FLOW 2-way 100% silicone catheter integral blue, central open tip, with plug <b>new</b> | page 29  |
| BLUECATH 2-way integral 100% silicone catheter, short, with plug                             | page 30  |
| CARE FLOW Nephrostomy catheter                                                               | page 31  |
| CARE FLOW 2-way integral 100% silicone catheter, with dilatation tip                         | page 32  |
|                                                                                              |          |
| <b>CARE FLOW Nephrostomy catheter</b>                                                        |          |
| CARE FLOW 2-way integral 100% silicone catheter, with dilatation tip                         | page 33  |

Table of Contents

|                                                                             |          |
|-----------------------------------------------------------------------------|----------|
| <b>CARE SET Catheter insertion kit and additional equipment</b>             |          |
| Catheter insertion kit Basic and Clinic                                     | page 34  |
| Catheter insertion kit Comfort and Comfort Plus                             | page 35  |
| Catheter insertion kit Comfort Eco                                          | page 36  |
| Additional equipment                                                        | page 37  |
|                                                                             |          |
| <b>CARE /ISC Intermittent Catheter</b>                                      |          |
| CARE CATH <i>Leir</i> IC disposable catheter, hydrophilic, 17 cm <b>new</b> | Seite 38 |
| CARE CATH <i>Leir</i> IC disposable catheter, hydrophilic, 41 cm <b>new</b> | Seite 49 |
|                                                                             |          |
| <b>CARE CON Male External Catheter</b>                                      |          |
| CARE CON Male External Catheter comfort / type I <b>new</b>                 | page 40  |
| CARE CON Male External Catheter special / type II <b>new</b>                | page 41  |
| CARE CON Male External Catheter active / type III <b>new</b>                | page 42  |
|                                                                             |          |
| <b>CARE FLOW Urine drainage systems and accessories</b>                     |          |
| CARE FLOW UB 4000 TUR Closed TUR irrigation bag system                      | page 43  |
| CARE FLOW UB 4000 Closed urine drainage system                              | page 44  |
| CARE FLOW UB 2000 Closed urine drainage system, incl. CARE BAG UB 2000      | page 45  |
| CARE FLOW UB 2000A Sterile single-use bedside bag                           | page 46  |
| CARE FLOW UB 2000E Sterile single-use bedside bag, incl. UB 2000US          | page 47  |
| CARE FLOW UB 2000USo, UB 1500USo, SB 1500USo                                | page 48  |
| CARE FLOW UB 1000 Closed urine drainage system for day & night              | page 49  |
| CARE FLOW UB 750/750US 1-chamber leg bag system                             | page 50  |
| CARE FLOW UB 750F/750UF 1-chamber leg bag systems with flip valve           | page 51  |
| CARE FLOW UB 500/500US 3-chamber leg bag system                             | page 52  |
|                                                                             |          |
| <b>CARE FIX Leg bag fixation</b>                                            |          |
| CARE FIX Universal Velcro strap                                             | page 53  |
| CARE FIX Velcro strap premium set, with and without buttons                 | page 54  |
| CARE FIX Leg pants with integrated leg bag pocket                           | page 55  |
| CARE FIX Fabric leg bag holder for the thigh                                | page 56  |
| CARE FIX Lower leg cuffs                                                    | page 57  |
| CARE FIX Premium one leg pants                                              | page 58  |
|                                                                             |          |
| <b>CARE FLOW Catheter valve</b>                                             |          |
| CARE FLOW and CARE VENT                                                     | page 59  |
|                                                                             |          |
| <b>CARE DRESS Dressing material</b>                                         |          |
| Suprapubic dressing set                                                     | page 60  |
| Abdominal PEG dressing set                                                  | page 61  |

German Health Care GmbH stands for innovative solutions in medical care of persons with physical disabilities. As an importer and distributor of products especially for the bladder management, we offer clinics, homecare providers and pharmacies a broad portfolio and individual solutions to improve the quality of life of your customers – whilst giving you the added safety of receiving high product quality as well as reliability.

Our own, extensive product range is thanks to the strategic cooperation with well-known manufacturers supplemented by established products. Due to our continuous exchange with our customers and the end-users, our products are not only in a state of steady improvement, but also new products are added to our assortment. Hence not only market-leading homecare providers value our high-quality product and services since several years, but also clinics and the purchasing communities of medical trade suppliers.

The production of GHC-products is guided by clear specifications of our manufacturers, based at different locations worldwide. To ensure the quality of the production processes, every production site is getting audited regularly. The manufacturers conduct quality controls directly after the production and before the goods gets shipped. The final inspections are conducted when the goods arrive at our high-bay warehouse in Wees. Before the GHC products are shipped to our customers, they get checked and tested extensively. Only perfectly produced products are getting released and our customers can rely on this strict quality policy.

This strict quality policy is also embedded in our quality management system according to DIN EN ISO 13485.

Not only can you rely on the quality of our products, but rather also on our employees, as they are the basis of our success. Our qualified medical device consultants in the back office as well as our field staff will gladly advise you on all questions concerning the GHC product portfolio and take care of it to process your requests as quickly as possible. Through continuous in-house workshops and further education as well as IT-supported processes we are guaranteeing this expertise, so that your inquiries, your orders, and your goods are processed as quickly as possible. You can count on that.

Ensure your incontinence supply quality by cooperating with GHC German Health Care GmbH, your system provider for revulsive incontinence care.



**Size & Color Coding: 1 Charrière (CH = 1 French 0.33mm)**

We carry catheters of different materials with different specifications to meet all the needs of our customers. Among them are catheters made of full silicone and Silicone elastomer for long-term drainage, but also latex catheter for short-term drainage.

## Classic catheter materials



100% Silicone



*Silicone elastomer –  
Latex core with  
silicone coating*



*Latex*

## Catheter tips for transurethral and suprapubic urine drainage



*Nelaton tip*



*Tiemann tip*



*central open tip*



*dilatation tip*



*Couvelaire tip*



*Mercier tip*



*Dufour tip*



*Flute tip*

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone catheter  
2-way balloon catheter with nelaton tip

Product specifications:

- ✓ In CH sizes 06 – 10 the catheter comes with a stabilising guiderail.
- ✓ Due to the high balloon capacity of up to 30 ml, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for long-term drainage.



Nelaton tip



Balloon volume: 03-30 ml

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 06 | 03 ml   | 27 cm  | 01/01/11/06/03/GHC | 01150196 |
| CH 08 | 05 ml   | 27 cm  | 01/01/11/08/05/GHC | 01150204 |
| CH 10 | 05 ml   | 27 cm  | 01/01/11/10/05/GHC | 01150210 |
| CH 12 | 10 ml   | 40 cm  | 01/01/11/12/10/GHC | 01150368 |
| CH 14 | 10 ml   | 40 cm  | 01/01/11/14/10/GHC | 01150374 |
| CH 16 | 30 ml   | 40 cm  | 01/01/11/16/30/GHC | 01150380 |
| CH 18 | 30 ml   | 40 cm  | 01/01/11/18/30/GHC | 01150397 |
| CH 20 | 30 ml   | 40 cm  | 01/01/11/20/30/GHC | 01150428 |
| CH 22 | 30 ml   | 40 cm  | 01/01/11/22/30/GHC | 01150434 |
| CH 24 | 30 ml   | 40 cm  | 01/01/11/24/30/GHC | 01150440 |
| CH 26 | 30 ml   | 40 cm  | 01/01/11/26/30/GHC | 01150457 |

Packaging Unit: 10 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW  
Transurethral *Soft* 100% silicone catheter  
2-way balloon catheter with nelaton tip

Product specifications:

- ✓ In CH sizes 06 – 10 the catheter comes with a stabilising guiderail.
- ✓ Very high wearing comfort, due to the soft-silicone material of the catheter shaft.
- ✓ Safe fixation of the catheter in the bladder due to the balloon.
- ✓ Use our blocking syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ The catheter can be used for long-term drainage.



Nelaton tip



Balloon volume: 03-30 ml

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 06 | 03 ml   | 31 cm  | 01/01/01/06/03/GHC | 17883185 |
| CH 08 | 05 ml   | 31 cm  | 01/01/01/08/05/GHC | 17883216 |
| CH 10 | 05 ml   | 31 cm  | 01/01/01/10/05/GHC | 17883222 |
| CH 12 | 10 ml   | 40 cm  | 01/01/01/12/10/GHC | 17883239 |
| CH 14 | 10 ml   | 40 cm  | 01/01/01/14/10/GHC | 17883251 |
| CH 16 | 30 ml   | 40 cm  | 01/01/01/16/30/GHC | 17883268 |
| CH 18 | 30 ml   | 40 cm  | 01/01/01/18/30/GHC | 17883274 |
| CH 20 | 30 ml   | 40 cm  | 01/01/01/20/30/GHC | 17883280 |
| CH 22 | 30 ml   | 40 cm  | 01/01/01/22/30/GHC | 17883297 |
| CH 24 | 30 ml   | 40 cm  | 01/01/01/24/30/GHC | 17883305 |
| CH 26 | 30 ml   | 40 cm  | 01/01/01/26/30/GHC | 17883311 |

Packaging Unit: 10 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW

Transurethral *integral* 100% silicone catheter

2-way integral balloon catheter with nelaton tip and plug

- Product specifications:
- ✓ Safe fixation of the catheter, due to the seamlessly incorporated *integral* balloon into the shaft.
  - ✓ Use our blocker syringe containing 10% glycerine solution here.
  - ✓ Excellent biocompatibility and high resistance to encrustation.
  - ✓ The X-ray line allows for radiological imaging of the catheter.
  - ✓ MRT – compatibility.
  - ✓ The catheter can be used for long-term drainage.



Nelaton tip,  
*integral* balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 41 cm  | 01/01/25/12/05/GHC | 18064437 |
| CH 14 | 05 ml   | 41 cm  | 01/01/25/14/05/GHC | 18064443 |
| CH 16 | 10 ml   | 41 cm  | 01/01/25/16/10/GHC | 18064466 |
| CH 18 | 10 ml   | 41 cm  | 01/01/25/18/10/GHC | 18064472 |
| CH 20 | 10 ml   | 41 cm  | 01/01/25/20/10/GHC | 18064489 |
| CH 22 | 10 ml   | 41 cm  | 01/01/25/22/10/GHC | 18064495 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml


Pictures may differ from the original

BLUECATH

Transurethral *integral* 100% silicone catheter

2-way integral balloon catheter with nelaton tip and plug

- Product specifications:
- ✓ Safe fixation of the catheter, due to the seamlessly incorporated *integral* balloon into the shaft.
  - ✓ Use our blocker syringe containing 10% glycerine solution here.
  - ✓ Excellent biocompatibility and high resistance to encrustation.
  - ✓ MRT – compatibility.
  - ✓ The catheter can be used for long-term drainage.



Nelaton tip,  
*integral* balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 41 cm  | 01/01/26/12/05/GHC | 18064503 |
| CH 14 | 05 ml   | 41 cm  | 01/01/26/14/05/GHC | 18064526 |
| CH 16 | 10 ml   | 41 cm  | 01/01/26/16/10/GHC | 18064532 |
| CH 18 | 10 ml   | 41 cm  | 01/01/26/18/10/GHC | 18064549 |
| CH 20 | 10 ml   | 41 cm  | 01/01/26/20/10/GHC | 18064555 |
| CH 22 | 10 ml   | 41 cm  | 01/01/26/22/10/GHC | 18064561 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original

CARE FLOW  
Transurethral silicone elastomer catheter  
2-way balloon catheter with nelaton tip

Product specifications:

- ✓ High elasticity and wearing comfort, due to the silicone coating on the inner and outer catheter wall.
- ✓ Due to the high balloon capacity, the catheter is safety fixed in the bladder.
- ✓ For blocking use our blocker syringe containing 10% glycerine solution.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The catheter can be used for long-term drainage.



Nelaton tip



Silicone elastomer –  
Latex core with silicone coating



Balloon volume: 10-30 ml

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 10 ml   | 40 cm  | 01/01/10/12/10/GHC | 01149371 |
| CH 14 | 30 ml   | 40 cm  | 01/01/10/14/30/GHC | 01149388 |
| CH 16 | 30 ml   | 40 cm  | 01/01/10/16/30/GHC | 01149394 |
| CH 18 | 30 ml   | 40 cm  | 01/01/10/18/30/GHC | 01149939 |
| CH 20 | 30 ml   | 40 cm  | 01/01/10/20/30/GHC | 01149945 |
| CH 22 | 30 ml   | 40 cm  | 01/01/10/22/30/GHC | 01149951 |
| CH 24 | 30 ml   | 40 cm  | 01/01/10/24/30/GHC | 01149968 |
| CH 26 | 30 ml   | 40 cm  | 01/01/10/26/30/GHC | 01150173 |

Packaging Unit: 10 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone catheter  
2-way balloon catheter with tiemann tip

Product specifications:

- ✓ An atraumatic dilating ball tip ensures a smooth passage of the prostatic fossa.
- ✓ Safe fixation of the catheter due to its high balloon capacity of up to 30 ml.
- ✓ Use our blocking syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility
- ✓ The catheter can be used for long term drainage.



Tiemann tip



Balloon volume: 5-30 ml

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 40 cm  | 01/01/08/12/05/GHC | 17933670 |
| CH 14 | 20 ml   | 40 cm  | 01/01/08/14/20/GHC | 17933687 |
| CH 16 | 30 ml   | 40 cm  | 01/01/08/16/30/GHC | 17933693 |
| CH 18 | 30 ml   | 40 cm  | 01/01/08/18/30/GHC | 17933701 |
| CH 20 | 30 ml   | 40 cm  | 01/01/08/20/30/GHC | 17933718 |
| CH 22 | 30 ml   | 40 cm  | 01/01/08/22/30/GHC | 17933724 |
| CH 24 | 30 ml   | 40 cm  | 01/01/08/24/30/GHC | 17933730 |

Packaging Unit: 10 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone catheter  
2-way balloon catheter with flute tip

Product specifications:

- ✓ Excellent drainage performance, due to the central open tip, drainage eye, and high lumen, even for problematic urine findings.
- ✓ An atraumatic flute tip ensures an easy placement of the catheter.
- ✓ Safe fixation of the catheter due to its high balloon capacity of up to 30 ml.
- ✓ Use our blocking syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The X-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for long term drainage.



Central open flute tip

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 10 ml   | 40 cm  | 01/01/07/12/10/GHC | 17933747 |
| CH 14 | 10 ml   | 40 cm  | 01/01/07/14/10/GHC | 17933753 |
| CH 16 | 30 ml   | 40 cm  | 01/01/07/16/30/GHC | 17933776 |
| CH 18 | 30 ml   | 40 cm  | 01/01/07/18/30/GHC | 17933782 |
| CH 20 | 30 ml   | 40 cm  | 01/01/07/20/30/GHC | 17933799 |
| CH 22 | 30 ml   | 40 cm  | 01/01/07/22/30/GHC | 17933807 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 10-30 ml

Pictures may differ from the original

CARE FLOW  
Transurethral silicone elastomer catheter  
2-way balloon catheter with couvelaire tip

Product specifications:

- ✓ High elasticity and wearing comfort, due to the silicone coating on the inner and outer catheter wall.
- ✓ Optimal drainage and irrigation due to two staggered drainage eyes and open couvelaire tip.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Safe fixation of the catheter, due to the high balloon capacity of up to 60 ml.
- ✓ For blocking use our blocker syringe containing 10% glycerine solution.
- ✓ The catheter can be used for long-term drainage.



Couvelaire tip

Silicone elastomer – Latex core with silicone coating

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 18 | 60 ml   | 40 cm  | 01/01/15/18/60/GHC | 06434975 |
| CH 20 | 60 ml   | 40 cm  | 01/01/15/20/60/GHC | 06434952 |
| CH 22 | 60 ml   | 40 cm  | 01/01/15/22/60/GHC | 06434969 |
| CH 24 | 60 ml   | 40 cm  | 01/01/15/24/60/GHC | 09710080 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: up to 60 ml

Pictures may differ from the original

CARE FLOW  
Transurethral siliconised latex catheter  
2-way balloon catheter with nelaton tip

Product specifications:

- ✓ Made of high-quality, siliconized medical latex.
- ✓ The catheter offers high patient comfort.
- ✓ Due to the high balloon capacity, the catheter is safety fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The catheter is suited for short-term drainage of up to 5 days only.
- ✓ This product contains natural latex which may cause allergic reactions.



Nelaton tip

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 10 ml   | 40 cm  | 01/01/13/12/10/GHC | 01150701 |
| CH 14 | 10 ml   | 40 cm  | 01/01/13/14/10/GHC | 01150859 |
| CH 16 | 30 ml   | 40 cm  | 01/01/13/16/30/GHC | 01150865 |
| CH 18 | 30 ml   | 40 cm  | 01/01/13/18/30/GHC | 01150871 |
| CH 20 | 30 ml   | 40 cm  | 01/01/13/20/30/GHC | 01150888 |
| CH 22 | 30 ml   | 40 cm  | 01/01/13/22/30/GHC | 01150919 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 10-30 ml

Pictures may differ from the original

CARE FLOW  
Transurethral PVC single-use female catheter  
Transparent with funnel and nelaton tip, DEHP-free

Product specifications:

- ✓ Disposable female catheter with two staggered atraumatic drainage eyes.
- ✓ Rounded and closed tip.
- ✓ Soft PVC material.
- ✓ DEHP-free.
- ✓ The catheter is for single use only.



Nelaton tip

| Size  | Length | Article number     | PZN      |
|-------|--------|--------------------|----------|
| CH 08 | 20 cm  | 01/01/35/08/00/GHC | 18320975 |
| CH 10 | 20 cm  | 01/01/35/10/00/GHC | 18320981 |
| CH 12 | 20 cm  | 01/01/35/12/00/GHC | 18320998 |
| CH 14 | 20 cm  | 01/01/35/14/00/GHC | 18321006 |
| CH 16 | 20 cm  | 01/01/35/16/00/GHC | 18321035 |

Packaging Unit: 100 pieces  
Packaging: sterile



Laterally displaced eyes

Pictures may differ from the original

CARE FLOW

Transurethral PVC single-use male catheter

Transparent with funnel, nelaton and tiemann tip

Product specifications:

- ✓ Disposable male catheter with two staggered atraumatic drainage eyes.
- ✓ Rounded and closed Nelaton or Tiemann tip, for a smooth passage of the prostatic fossa.
- ✓ Soft PVC material.
- ✓ The catheter is for single use only.



Nelaton tip

Tiemann tip

Nelaton tip

| Size  | Length | Article number     | PZN      |
|-------|--------|--------------------|----------|
| CH 10 | 40 cm  | 01/01/20/10/00/GHC | 08871349 |
| CH 12 | 40 cm  | 01/01/20/12/00/GHC | 08871355 |
| CH 14 | 40 cm  | 01/01/20/14/00/GHC | 08871361 |
| CH 16 | 40 cm  | 01/01/20/16/00/GHC | 08871378 |
| CH 18 | 40 cm  | 01/01/20/18/00/GHC | 08871384 |

Tiemann tip

| Size  | Length | Article number     | PZN      |
|-------|--------|--------------------|----------|
| CH 12 | 40 cm  | 01/01/21/12/00/GHC | 01929620 |
| CH 14 | 40 cm  | 01/01/21/14/00/GHC | 01929637 |
| CH 16 | 40 cm  | 01/01/21/16/00/GHC | 01929666 |
| CH 18 | 40 cm  | 01/01/21/18/00/GHC | 01971792 |
| CH 20 | 40 cm  | 01/01/21/20/00/GHC | 09219740 |

Packaging Unit: CH 08-18: 50 pieces each, CH 20: 30 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW

Transurethral 100% silicone irrigation catheter

3-way balloon catheter with nelaton tip

Product specifications:

- ✓ An atraumatic Nelaton tip ensures a smooth passage of the prostatic fossa.
- ✓ Optimal drainage and irrigation performance, due to two staggered drainage eyes.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for post-operative long-term drainage.



Nelaton tip

Balloon volume: 30 ml

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 18 | 30 ml   | 40 cm  | 01/01/22/18/30/GHC | 15618237 |
| CH 20 | 30 ml   | 40 cm  | 01/01/22/20/30/GHC | 15618243 |
| CH 22 | 30 ml   | 40 cm  | 01/01/22/22/30/GHC | 15618266 |
| CH 24 | 30 ml   | 40 cm  | 01/01/22/24/30/GHC | 15618295 |

Packaging Unit: 10 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone irrigation catheter  
3-way balloon catheter with dufour tip

Product specifications:

- ✓ An atraumatic and kink resistant Defour tip ensures a smooth passage of the prostatic fossa.
- ✓ Optimal drainage and irrigation performance, due to two staggered drainage eyes and an open tip.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray tip allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for post-operative long-term drainage.



Dufour tip

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 18 | 40 ml   | 40 cm  | 01/01/16/18/40/GHC | 01151178 |
| CH 20 | 50 ml   | 40 cm  | 01/01/16/20/50/GHC | 01151184 |
| CH 22 | 60 ml   | 40 cm  | 01/01/16/22/60/GHC | 01151190 |
| CH 24 | 70 ml   | 40 cm  | 01/01/16/24/70/GHC | 01151379 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 40-70 ml

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone irrigation catheter  
3-way balloon catheter with couvelaire tip

Product specifications:

- ✓ Optimal drainage and irrigation due to two staggered drainage eyes and an open tip.
- ✓ An atraumatic and kink resistant Couvelaire tip ensures a smooth passage of the prostatic fossa.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray tip allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ The catheter can be used for post-operative long-term drainage.



Couvelaire tip

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 18 | 40 ml   | 40 cm  | 01/01/17/18/40/GHC | 01151385 |
| CH 20 | 50 ml   | 40 cm  | 01/01/17/20/50/GHC | 01151391 |
| CH 22 | 60 ml   | 40 cm  | 01/01/17/22/60/GHC | 01151416 |
| CH 24 | 70 ml   | 40 cm  | 01/01/17/24/70/GHC | 01151439 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 40-70 ml

Pictures may differ from the original

CARE FLOW  
Transurethral 100% silicone irrigation catheter  
3-way balloon catheter with mercier tip

- Product specifications:
- ✓ An atraumatic and kink resistant Mercier tip ensures a smooth passage of the prostatic fossa.
  - ✓ Optimal drainage and irrigation due to two staggered drainage eyes and an open tip.
  - ✓ Excellent biocompatibility and high resistance to encrustation.
  - ✓ Due to the high balloon capacity, the catheter is safely fixed in the bladder.
  - ✓ Use our blocker syringe containing 10% glycerine solution here.
  - ✓ The x-ray tip allows for radiological imaging of the catheter.
  - ✓ MRT – compatibility.
  - ✓ The catheter can be used for post-operative long-term drainage.



Mercier tip

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 18 | 40 ml   | 40 cm  | 01/01/18/18/40/GHC | 01151511 |
| CH 20 | 50 ml   | 40 cm  | 01/01/18/20/50/GHC | 01151528 |
| CH 22 | 60 ml   | 40 cm  | 01/01/18/22/60/GHC | 01151534 |
| CH 24 | 70 ml   | 40 cm  | 01/01/18/24/70/GHC | 01151540 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 40-70 ml

CARE FLOW UNIBAL  
Suprapubic *integral* 100% silicone catheter  
2-way integral balloon catheter

- Product specifications:
- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
  - ✓ The catheter offers high patient comfort.
  - ✓ High drainage performance, due to the high lumen and the central open tip.
  - ✓ We recommend using or blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter
  - ✓ Excellent biocompatibility and high resistance to encrustation.
  - ✓ Short and central open tip.
  - ✓ The x-ray line allows for radiological imaging of the catheter.
  - ✓ MRT – compatibility.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml   | 40 cm  | 01/02/12/10/03/GHC | 01152568 |
| CH 12 | 05 ml   | 40 cm  | 01/02/12/12/05/GHC | 01152574 |
| CH 14 | 05 ml   | 40 cm  | 01/02/12/14/05/GHC | 01152605 |
| CH 16 | 10 ml   | 40 cm  | 01/02/12/16/10/GHC | 01152611 |
| CH 18 | 10 ml   | 40 cm  | 01/02/12/18/10/GHC | 01152628 |
| CH 20 | 10 ml   | 40 cm  | 01/02/12/20/10/GHC | 01152700 |
| CH 22 | 10 ml   | 40 cm  | 01/02/12/22/10/GHC | 01152717 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 03-10 ml

Pictures may differ from the original

Pictures may differ from the original

CARE FLOW UNIBAL  
Suprapubic *integral* 100% silicone catheter  
2-way integral balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers high patient comfort.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ We recommend using or blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml   | 40 cm  | 01/02/14/10/03/GHC | 10198115 |
| CH 12 | 05 ml   | 40 cm  | 01/02/14/12/05/GHC | 10198121 |
| CH 14 | 05 ml   | 40 cm  | 01/02/14/14/05/GHC | 10198138 |
| CH 16 | 10 ml   | 40 cm  | 01/02/14/16/10/GHC | 10198144 |
| CH 18 | 10 ml   | 40 cm  | 01/02/14/18/10/GHC | 10198150 |
| CH 20 | 10 ml   | 40 cm  | 01/02/14/20/10/GHC | 10198167 |
| CH 22 | 10 ml   | 40 cm  | 01/02/14/22/10/GHC | 10198173 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 03-10 ml

Pictures may differ from the original

CARE FLOW  
Suprapubic 100% silicone catheter *integral*  
2-way balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers an excellent wearing comfort.
- ✓ We recommend using our blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter
- ✓ 100% silicone for excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 41 cm  | 01/02/09/12/05/GHC | 16840937 |
| CH 14 | 05 ml   | 41 cm  | 01/02/09/14/05/GHC | 16840943 |
| CH 16 | 10 ml   | 41 cm  | 01/02/09/16/10/GHC | 16840972 |
| CH 18 | 10 ml   | 41 cm  | 01/02/09/18/10/GHC | 16840989 |
| CH 20 | 10 ml   | 41 cm  | 01/02/09/20/10/GHC | 16840995 |
| CH 22 | 10 ml   | 41 cm  | 01/02/09/22/10/GHC | 16841003 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original

CARE FLOW UNIBAL  
Suprapubic *integral* 100% silicone catheter  
2-way integral balloon catheter, short

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml   | 33 cm  | 01/02/13/10/03/GHC | 01152746 |
| CH 12 | 05 ml   | 33 cm  | 01/02/13/12/05/GHC | 01152812 |
| CH 14 | 05 ml   | 33 cm  | 01/02/13/14/05/GHC | 01152829 |
| CH 16 | 10 ml   | 33 cm  | 01/02/13/16/10/GHC | 01152841 |
| CH 18 | 10 ml   | 33 cm  | 01/02/13/18/10/GHC | 01152918 |
| CH 20 | 10 ml   | 33 cm  | 01/02/13/20/10/GHC | 01152924 |
| CH 22 | 10 ml   | 33 cm  | 01/02/13/22/10/GHC | 01152947 |

Packaging Unit: 10 pieces  
Packaging: sterile




Balloon volume: 03-10 ml

Pictures may differ from the original

BLUECATH  
Suprapubic *integral* 100% silicone catheter  
2-way integral balloon catheter with dilatation tip and plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Dilatation tip with central bore.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optional placement.



Dilatation tip with central bore, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 40 cm  | 01/02/21/12/05/GHC | 11850948 |
| CH 14 | 05 ml   | 40 cm  | 01/02/21/14/05/GHC | 11850954 |
| CH 16 | 10 ml   | 40 cm  | 01/02/21/16/10/GHC | 11850983 |
| CH 18 | 10 ml   | 40 cm  | 01/02/21/18/10/GHC | 11851008 |
| CH 20 | 10 ml   | 40 cm  | 01/02/21/20/10/GHC | 11851014 |
| CH 22 | 10 ml   | 40 cm  | 01/02/21/22/10/GHC | 11851020 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original

BLUECATH  
Suprapubic *integral* 100% silicone catheter  
2-way integral balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ Short, central open tip.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optional placement.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 40 cm  | 01/02/22/12/05/GHC | 11851037 |
| CH 14 | 05 ml   | 40 cm  | 01/02/22/14/05/GHC | 11851043 |
| CH 16 | 10 ml   | 40 cm  | 01/02/22/16/10/GHC | 11851072 |
| CH 18 | 10 ml   | 40 cm  | 01/02/22/18/10/GHC | 11851089 |
| CH 20 | 10 ml   | 40 cm  | 01/02/22/20/10/GHC | 11851095 |
| CH 22 | 10 ml   | 40 cm  | 01/02/22/22/10/GHC | 11851126 |
| CH 24 | 10 ml   | 40 cm  | 01/02/22/24/10/GHC | 11851132 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original

CARE FLOW  
Suprapubic 100% silicone catheter *integral blue*  
2-way balloon catheter with plug

Product specifications:

- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
- ✓ The catheter offers an excellent wearing comfort.
- ✓ We recommend using our blocking syringe containing a 10% glycerine solution, for a safe fixation of the catheter.
- ✓ 100% silicone for excellent biocompatibility and high resistance to encrustation.
- ✓ Short and central open tip.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optimal placement.



Central open, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 41 cm  | 01/02/11/12/05/GHC | 16841026 |
| CH 14 | 05 ml   | 41 cm  | 01/02/11/14/05/GHC | 16841032 |
| CH 16 | 10 ml   | 41 cm  | 01/02/11/16/10/GHC | 16841049 |
| CH 18 | 10 ml   | 41 cm  | 01/02/11/18/10/GHC | 16841061 |
| CH 20 | 10 ml   | 41 cm  | 01/02/11/20/10/GHC | 16841078 |
| CH 22 | 10 ml   | 41 cm  | 01/02/11/22/10/GHC | 16841084 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

Pictures may differ from the original

BLUECATH

Suprapubic *integral* 100% silicone catheter

2-way integral balloon catheter with plug, short

- Product specifications:
- ✓ The catheter is equipped with an integral balloon, resulting in an even thickness throughout the catheter. This ensures that there is no painful compression of the fascia during catheter placement.
  - ✓ High drainage performance, due to the high lumen and the central open tip.
  - ✓ Safe fixation of the catheter, due to the *integral* balloon.
  - ✓ Use our blocker syringe containing 10% glycerine solution here.
  - ✓ Excellent biocompatibility and high resistance to encrustation.
  - ✓ Short, central open tip.
  - ✓ MRT – compatibility.
  - ✓ Graduation marks to ensure an optional placement.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 12 | 05 ml   | 33 cm  | 01/02/23/12/05/GHC | 11851149 |
| CH 14 | 05 ml   | 33 cm  | 01/02/23/14/05/GHC | 11851155 |
| CH 16 | 10 ml   | 33 cm  | 01/02/23/16/10/GHC | 11851161 |
| CH 18 | 10 ml   | 33 cm  | 01/02/23/18/10/GHC | 11851178 |
| CH 20 | 10 ml   | 33 cm  | 01/02/23/20/10/GHC | 11851184 |
| CH 22 | 10 ml   | 33 cm  | 01/02/23/22/10/GHC | 11851190 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 05-10 ml

CARE FLOW UNIBAL

Suprapubic puncture set

2-way *integral* balloon catheter

- Product specifications:
- The puncture set contains:
- 100% silicone integral balloon catheter
  - Stainless steel cannula with special cut
  - Scalpel with small handle
- ✓ *Integral* balloon catheter with seamlessly incorporated balloon into the shaft.
  - ✓ Optimal use of the inside diameter of the cannula due to integral balloon catheter: Charrière 10 and 12 with split cannula CH 13.8 (4.6 mm), puncture length 120 mm. Charrière 14 with split cannula CH 16.8 (5.6 mm), puncture length 120 mm.
  - ✓ High drainage performance, due to the high lumen and the central open tip.
  - ✓ Safe fixation of the catheter, due to the integral balloon.
  - ✓ Use our blocker syringe containing 10% glycerine solution here.
  - ✓ The x-ray line allows for radiological imaging of the catheter.
  - ✓ MRT – compatibility.



Central open tip, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml   | 40 cm  | 01/02/15/10/03/GHC | 01153059 |
| CH 12 | 05 ml   | 40 cm  | 01/02/15/12/05/GHC | 01153065 |
| CH 14 | 05 ml   | 40 cm  | 01/02/15/14/05/GHC | 01153088 |

Packaging Unit: 5 pieces  
Packaging: sterile



Balloon volume: 03-05 ml

Pictures may differ from the original

Pictures may differ from the original

CARE FLOW UNIBAL  
Suprapubic change set, incl. guidewire  
2-way integral balloon catheter

Product specifications:

This charge set contains:

- 100% silicone *integral* balloon catheter
- Teflon-coated guidewire

- ✓ Integral balloon catheter with seamlessly incorporated balloon into the shaft.
- ✓ High drainage performance, due to the high lumen and the central open tip.
- ✓ Stable teflon-coated guidewire.
- ✓ Safe fixation of the catheter, due to the integral balloon.
- ✓ Use our blocker syringe containing 10% glycerine solution here.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ Teflon-coated guidewire, lenght 80 cm, 3 cm straight and flexible tip, thickness 0,035”.



Central open tip, integral balloon

Guidewire

| Wire   | Tip        | Length | Article number     | PZN      |
|--------|------------|--------|--------------------|----------|
| coated | 3 cm flex. | 80 cm  | 01/02/18/80/35/GHC | 14166431 |

Packaging Unit: 1 piece  
Packaging: sterile

Change set, incl. Guidewire

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 03 ml   | 40 cm  | 01/02/16/10/03/GHC | 01153154 |
| CH 12 | 05 ml   | 40 cm  | 01/02/16/12/05/GHC | 01153160 |
| CH 14 | 05 ml   | 40 cm  | 01/02/16/14/05/GHC | 01153183 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 03-05 ml

NEPHROSTOMY  
Nephrostomy *integral* 100% silicone catheter  
2-way integral balloon catheter with dilatation tip

Product specifications:

- ✓ Excellent biocompatibility and high resistance to encrustation.
- ✓ The catheters are suitable for long-term drainage.
- ✓ Dilating nelaton tip with central bore.
- ✓ High drainage performance, due to four drainage eyes located in front of and behind the balloon.
- ✓ The x-ray line allows for radiological imaging of the catheter.
- ✓ MRT – compatibility.
- ✓ Graduation marks to ensure an optimal placement.
- ✓ 100% silicone.



Dilating nelaton tip with central bore, integral balloon

| Size  | Balloon | Length | Article number     | PZN      |
|-------|---------|--------|--------------------|----------|
| CH 10 | 1,5 ml  | 40 cm  | 01/03/10/10/01/GHC | 13851171 |
| CH 12 | 2,0 ml  | 40 cm  | 01/03/10/12/02/GHC | 13851188 |
| CH 14 | 3,0 ml  | 40 cm  | 01/03/10/14/03/GHC | 13851202 |
| CH 16 | 5,0 ml  | 40 cm  | 01/03/10/16/05/GHC | 13851219 |
| CH 18 | 5,0 ml  | 40 cm  | 01/03/10/18/05/GHC | 13851225 |

Packaging Unit: 10 pieces  
Packaging: sterile



Balloon volume: 1,5 - 5,0 ml

Pictures may differ from the original

Pictures may differ from the original

CARE SET  
Catheter insertion kit Basic and Clinic

Product specifications:

- ✓ The Basic Set comprises an inner and an outer package.
- ✓ The outer package of the basic set contains:  
one sterile 10 ml syringe for unblocking the placed catheter.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



| Article | Article number   | PZN      |
|---------|------------------|----------|
| Basic   | 01/01/99/BS/GHC  | 01972188 |
| Clinic  | 01/01/99/BSK/GHC | 08845872 |

CARE SET  
Catheter insertion kit Comfort and Comfort Plus

Product specifications:

- ✓ The Catheter insertion Kits comprises an inner and an outer package.
- ✓ The outer package of the basic set contains: one sterile 10 ml syringe for unblocking the inserted catheter a mucous membrane disinfectant Octenisept 15 ml.
- ✓ Additional included in Comfort Plus: lubricant Endosgel 6 ml.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



| Article      | Article number     | PZN      |
|--------------|--------------------|----------|
| Comfort      | 01/01/99/KL/CO/GHC | 14209815 |
| Comfort Plus | 01/01/99/KL/CP/GHC | 14209821 |

Components of the catheter insertion kit are:

1 Tray



6 Swabs



1 Pair of tweezers



2 Underpads



1 Blocker syringe with 10% glycerine



3 Nitrile gloves



1 Fenestrated drape



Packaging Unit: 20 pieces  
Packaging: sterile

Additional included in Basic set:



10 ml syringe

Pictures may differ from the original

Components of the catheter insertion kit are:


1 Tray



6 Swabs



1 Pair of tweezers



10 ml Syringe




1 Blocker syringe with 10% glycerine



3 Nitrile gloves



1 Fenestrated drape



2 Underpads



Mucous membrane disinfectant



Octenisept 15 ml

Additional included in Comfort Plus:



Lubricant Endosgel 6 ml

Packaging Unit: 35 pieces  
Packaging: sterile

Pictures may differ from the original

CARE SET  
Catheter insertion kit Comfort Eco

Product specifications:

- ✓ The Comfort Eco Set comprises an inner and an outer package.
- ✓ The outer package of the basic set contains: one sterile 10 ml syringe for unblocking the inserted catheter, mucous membrane disinfectant Octenisept 15 ml, lubricant LUBRICARE 6 ml.
- ✓ The items in the inner, sterilised package are used to place a new transurethral catheter.



| Article     | Article number     | PZN      |
|-------------|--------------------|----------|
| Comfort Eco | 01/01/99/KL/CE/GHC | 16134211 |

Additional equipment

Accessories for the catheter replacement

Next to our assortment covering catheters as well as drainage- and holding- systems, we also provide necessary accessories, to maintain a satisfactory supply of revulsive incontinence aids and ensuring an optimal replacement of a catheter, which all helps to maintain a high level of supply quality.

We want to highlight:

- The **blocking syringe containing a 10% glycerine solution, 10 ml** (Art.-Nr. 01/01/99/KB/10/GHC), ensures a steady volume in the catheter balloon during the wearing period, by preventing the liquid medium from diffusing out of the balloon. Thereby the balloon catheter stays fixated, and the risk of a dislocated catheter can almost be excluded.
- The **lubricant LUBRICARE, 6 ml** (Art.-Nr. 01/01/99/GG/001), adhere excellently onto the mucous membrane and has excellent lubrication properties, which ensures a smooth replacement of a catheter.

The below illustrated products are only an extract of our comprehensive assortment. We gladly provide more information upon request.




Blocker syringe with 10% glycerine, 10 ml



Lubricant LUBRICARE, 06 ml



Catheter plug



Octenisept, 15 ml

Components of the catheter insertion kit are:

|                                      |                  |                     |               |                   |
|--------------------------------------|------------------|---------------------|---------------|-------------------|
| 1 Tray                               | 6 Swabs          | 1 Pair of tweezers  | 10 ml Syringe | Mucous membrane   |
|                                      |                  |                     |               |                   |
| 1 Blocker syringe with 10% glycerine | 3 Nitrile gloves | 1 Fenestrated drape | 2 Underpads   | Octenisept, 15 ml |
|                                      |                  |                     |               |                   |
|                                      |                  |                     |               | Lubricant         |
|                                      |                  |                     |               |                   |
|                                      |                  |                     |               | LUBRICARE, 06 ml  |

Packaging Unit: 35 pieces  
Packaging: sterile

| Article                | Article number      | PZN      | VE                  |
|------------------------|---------------------|----------|---------------------|
| 1. Blocker syringe     | 01/01/99/KB/10/GHC  | 06718980 | 1VE = 200 pieces    |
| 2. LUBRICARE Lubricant | 01/01/99/GG/06/001  | 13360220 | 1VE = 100 pieces    |
| LUBRICARE Lubricant    | 01/01/99/GG/06/25   | 16938435 | 1VE = 4 x 25 pieces |
| Endosgel Lubricant     | 01/01/99/HG/EG/GHC  | 02437623 | 1VE = 250 pieces    |
| 3. Catheter plug       | 01/01/99/ZB/ST/GHC  | 14354919 | 1VE = 20 pieces     |
| 4. Octenisept          | 01/01/99/OCT/15/GHC | 04804008 | 1VE = 100 pieces    |

Packaging: sterile

Pictures may differ from the original

Pictures may differ from the original

CARE CATH<sub>Leiv</sub>

IC disposable catheter, hydrophilic, 17 cm

Ready for immediate use, Nelaton

Product specifications:

- ✓ Very good gliding properties thanks to hydrophilic coating.
- ✓ Ready for immediate use thanks to ready-to-use, hydrophilic coating.
- ✓ Atraumatic application due to rounded drainage eyes.
- ✓ Optimal length of 17 cm ensures complete bladder emptying.
- ✓ Non-contact application and good guideability due to the correct rigidity of the catheter.
- ✓ Good grip of the handle allows good handling, even with limited dexterity.
- ✓ Easy connection of urine collection bag and/or instillation solution via the handle cone.



Atraumatic, rounded drainage eyes



Grip of the handle

CARE CATH<sub>Leiv</sub>

| Size  | Length | Item      | number   |
|-------|--------|-----------|----------|
| CH 10 | 17 cm  | LIV-H1017 | 18890885 |
| CH 12 | 17 cm  | LIV-H1217 | 18890916 |
| CH 14 | 17 cm  | LIV-H1417 | 18890922 |

Packaging unit: 30 pieces  
Packaging: sterile

Illustrations may differ from the original

CARE CATH<sub>Leiv</sub>

IC disposable catheter, hydrophilic, 41 cm

Ready for immediate use, Nelaton

Product specifications:

- ✓ Very good gliding properties thanks to hydrophilic coating.
- ✓ Ready for immediate use thanks to ready-to-use, hydrophilic coating.
- ✓ Atraumatic application due to rounded drainage eyes.
- ✓ Non-contact application and good guideability due to the correct rigidity of the catheter.
- ✓ Good grip of the handle allows good handling, even with limited dexterity.
- ✓ Easy connection of urine collection bag and/or instillation solution via the handle cone.



Atraumatic, rounded drainage eyes



Grip of the handle

CARE CATH<sub>Leiv</sub>

| Size  | Length | Item       | number   |
|-------|--------|------------|----------|
| CH 12 | 41 cm  | LEIV-H1241 | 18890833 |
| CH 14 | 41 cm  | LEIV-H1441 | 18890856 |
| CH 16 | 41 cm  | LEIV-H1641 | 18890862 |

Packaging unit: 30 pieces  
Packaging: sterile

Pictures may differ from the original

CARE CON comfort | Typ I  
Male External Catheter (MEC) 100% silicone  
Self-adhesive

Product specifications:

- ✓ The CARE CON comfort captivates through its comfortable silicone material and good hold.
- ✓ The material thickness of the urinary sheath has been developed to cling softly and flexibly whilst ensuring a secure hold, which provides high safety.
- ✓ The stable connecting tube prevents the MEC from bending and ensures a secure connection to the urine bag.
- ✓ Small steps between the sizes ensure perfect fit.
- ✓ The medical grade silicone maintains the breathing activity of the skin.
- ✓ Due to the use of transparent materials, skin alterations or irritations can be observed at any time.
- ✓ With a shaft length of 8,0 cm and an adhesive area of 4,5 cm, the CARE CON comfort meets the requirements of most MEC-users.



CARE CON special | Type II  
Male External Catheter (MEC) 100% silicone  
Self-adhesive

Product specifications:

- ✓ The CARE CON special captivates through its comfortable silicone material and good hold, despite the short shaft length.
- ✓ The material thickness of the urinary sheath has been developed to cling softly and flexibly whilst ensuring a secure hold, which provides high safety.
- ✓ The stable connecting tube prevents the MEC from bending and ensures a secure connection to the urine bag.
- ✓ Small steps between the sizes ensure perfect fit.
- ✓ With a shaft length of 5,0 cm and an adhesive area of 4,0 cm, the CARE CON special meets the special requirements of users with a retracted penis.
- ✓ The medical grade silicone maintains the breathing activity of the skin.
- ✓ Due to the use of transparent materials, skin alterations or irritations can be observed at any time.



| Diameter | Shaft / Adhesive surface | Article number     | PZN      |
|----------|--------------------------|--------------------|----------|
| 25 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/25/GHC | 17520554 |
| 28 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/28/GHC | 17520583 |
| 30 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/30/GHC | 17520620 |
| 32 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/32/GHC | 17520637 |
| 36 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/36/GHC | 17520643 |
| 41 mm    | 8,0 cm / 4,5 cm          | 01/01/28/01/41/GHC | 17520666 |

Packaging Unit: 30 pieces  
Packaging: non-sterile

Pictures may differ from the original

| Diameter | Shaft / Adhesive surface | Article number     | PZN      |
|----------|--------------------------|--------------------|----------|
| 25 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/25/GHC | 17520672 |
| 28 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/28/GHC | 17520689 |
| 30 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/30/GHC | 17520703 |
| 32 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/32/GHC | 17520726 |
| 36 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/36/GHC | 17520732 |
| 41 mm    | 5,0 cm / 4,0 cm          | 01/01/28/02/41/GHC | 17520749 |

Packaging Unit: 30 pieces  
Packaging: non-sterile

Pictures may differ from the original

CARE CON active | Typ III

Male External Catheter (MEC) 100% silicone

Self-adhesive

Product specifications:

- ✓ The CARE CON active captivates through its comfortable silicone material and good hold, with an extra-long adhesive area.
- ✓ The material thickness of the urinary sheath has been developed to cling softly and flexibly whilst ensuring a secure hold, which provides high safety.
- ✓ The stable connecting tube prevents the MEC from bending and ensures a secure connection to the urine bag.
- ✓ Small steps between the sizes ensure perfect fit.
- ✓ With a shaft length of 8,0 cm and an adhesive area of 7,0 cm, the CARE CON active meets the requirements of sporty and active MEC-users.
- ✓ The medical grade silicone maintains the breathing activity of the skin.
- ✓ Due to the use of transparent materials, skin alterations or irritations can be observed at any time.



| Diameter | Shaft / Adhesive surface | Article number     | PZN      |
|----------|--------------------------|--------------------|----------|
| 25 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/25/GHC | 17520755 |
| 28 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/28/GHC | 17520761 |
| 30 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/30/GHC | 17520784 |
| 32 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/32/GHC | 17520790 |
| 36 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/36/GHC | 17520809 |
| 41 mm    | 8,0 cm / 7,0 cm          | 01/01/28/03/41/GHC | 17520815 |

Packaging Unit: 30 pieces

Packaging: non-sterile

Pictures may differ from the original

CARE FLOW UB 4000TUR closed TUR irrigation bag system

with pump ball for controlled post-operative bladder irrigation

Capacity: 4000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder with fixing cord.
- ✓ Bacteria- and fluid-tight venting filters ensure permanent pressure compensation and provide for optimal drainage in the system.
- ✓ Due to controlled irrigation of the bladder with the pump ball integrated in the feeding tube, post-operative tamponade of the bladder can be avoided.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



| Capacity | Tube length | Article number     | PZN      |
|----------|-------------|--------------------|----------|
| 4000 ml  | 120 cm      | 04/02/4000/TUR/GHC | 06189004 |



Pump ball

The specific interaction of the pump ball and the two tube clamps allows to flush out blood clots in the bladder and tube system.



Drip chamber

An anti-reflux valve divides the drip chamber and reliably prevents a backflow of urine and thereby obviates the risk against an ascending infection.



Profiled tube

The octagonal inner profile of the feeding tube assures excellent stability, but at the same time it is also flexible. A high lumen guarantees a safe drainage into the drip chamber in any situation.



Drain tap in the doghouse

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

Packaging Unit: 10 pieces

Packaging: sterile

Pictures may differ from the original

CARE FLOW UB 4000  
Closed urine drainage system  
Capacity: 4000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder with fixing cord.
- ✓ Bacteria- and fluid-tight venting filters ensure permanent pressure compensation and provide for optimal drainage in the system.
- ✓ Clear graduation markings allow for reliable volume control.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.
- ✓ T-drain tap.
- ✓ Kink resistant tube.



| Capacity | Tube length | Article number     | PZN      |
|----------|-------------|--------------------|----------|
| 4000 ml  | 120 cm      | 04/02/4000/120/GHC | 06718939 |

CARE FLOW UB 2000  
Closed urine drainage system  
Capacity: 2000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder with fixing cord.
- ✓ Bacteria- and fluid-tight venting filters ensure permanent pressure compensation and provide for optimal drainage in the system.
- ✓ Clear graduation markings allow for reliable volume control.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.
- ✓ T-drain tap.
- ✓ Kink resistant tube.



| Capacity | Tube length | Article number     | PZN      |
|----------|-------------|--------------------|----------|
| 2000 ml  | 120 cm      | 04/02/2000/120/GHC | 06718922 |
| 2000 ml  | 90 cm       | 04/02/2000/90/GHC  | 09151597 |

**Packaging Unit:** 20 pieces  
**Packaging:** sterile



**Staggered connector**

The staggered connector provides a reliable connection with the catheter. A needle-free sampling of urine is possible at any time, due to the needle-free sample port, which prevents needlestick injuries.



**Profiled tube**

The octagonal inner profile of the feeding tube assures excellent stability, but at the same time it is also flexible. A high lumen guarantees a safe drainage into the drip chamber in any situation.



**Drip chamber**

An anti-reflux valve divides the drip chamber and reliably prevents a backflow of urine and thereby obviates the risk against an ascending infection.



**Drain tap in the doghouse**

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

**Packaging Unit:** 20 pieces  
**Packaging:** sterile

Pictures may differ from the original

The **CARE FLOW BAG** is a cover for the **CARE FLOW UB 2000** serving as a discreet way to store the urinary bedside bag.

Product specifications:

- ✓ Thanks to the convenient Velcro fastening the cover can be affixed to the urinary bag quickly and easily. A second Velcro fastening on the front allows for a clean and uncomplicated extraction and reinsertion of the drain tap.
- ✓ The cover is made from 100% polypropylene (PP) and is thus resistant against organic acids.
- ✓ The cover is light grey, washable at 30° and hence reusable.



| Article number   | PZN      |
|------------------|----------|
| 04/99/01/945/GHC | 09764331 |

**Packaging Unit:** 10 pieces  
**Packaging:** non-sterile

Pictures may differ from the original

CARE FLOW UB 2000A Disposable bedside bag  
Sterile bedside bag with integrated twin holder  
Capacity: 2000 ml

Product specifications:

- ✓ Can be fixed safely and correctly to many different types of bed frames by means of the twin holder.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.
- ✓ Clear graduation markings allow for reliable volume control.
- ✓ A slit ensures the safe and hygienic positioning of the dRAINTap.
- ✓ T-drain tap.
- ✓ Kink resistant tube.



| Capacity | Tube length | Article number     | PZN      |
|----------|-------------|--------------------|----------|
| 2000 ml  | 110 cm      | 04/04/2000/110/GHC | 08871289 |

CARE FLOW UB 2000E, UB 2000US  
Sterile / non-sterile single-use bedside bag with holder  
Capacity: 2000 ml

Product specifications:

- ✓ The system is equipped with a universal staggered cone and can be connected to all common catheters.
- ✓ The product must be changed daily.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.
- ✓ Use the twin holder shown below to attach the bag to many different types of bed frames.
- ✓ The urine bag is emptied safely by means of the T-valve.



| Article   | Capacity | Sterility   | Tube length | Article number     | PZN      |
|-----------|----------|-------------|-------------|--------------------|----------|
| UB 2000E  | 2000 ml  | sterile     | 90 cm       | 04/04/2000/90/GHC  | 16047274 |
| UB 2000US | 2000 ml  | non-sterile | 90 cm       | 04/06/2000/90/GHC  | 12418549 |
| UB 2000US | 2000 ml  | non-sterile | 130 cm      | 04/06/2000/130/GHC | 08845487 |

UB 2000E: **Packaging Unit:** 100 pieces, **Packaging:** sterile  
UB 2000US: **Packaging Unit:** 100 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces

Holder for UB 2000E and UB 2000US

- ✓ Plastic hanger for bedside bags with twin holder. Equipped with two pins taking up single-use bedside bags hanger holes.

| Article number    | PZN      |
|-------------------|----------|
| 04/99/BHL/PVC/GHC | 06718951 |

**Packaging Unit:** 100 pieces **Sub-packaging:** 10 pieces each  
**Packaging:** non-sterile



46



Staggered connector

The staggered connector provides a reliable connection with the catheter. A needle-free sampling of urine is possible at any time, due to the needle-free sample port, which prevents needlestick injuries.



Profiled tube

The octagonal inner profile of the feeding tube assures excellent stability, but at the same time it is also flexible. A high lumen guarantees a safe drainage into the drip chamber in any situation.



Venting filter

A bacteria and liquid-tight venting filters ensure permanent pressure equalisation. This ensures an optimal drainage of the system.



Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

**Packaging Unit:** 50 pieces  
**Packaging:** sterile

Pictures may differ from the original

47

Pictures may differ from the original

CARE FLOW UB 2000USo, UB 1500USo, SB 1500USo  
Non-sterile single-use bedside bag/secretion bag without drain tap  
Capacity: 2000 ml, 1500 ml

Product specifications:

- ✓ The system is equipped with a universal staggered cone and can be connected to all common catheters.
- ✓ The product must be changed daily.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria. (except the SB 1500USo)
- ✓ Use the twin holder shown below to attach the bag to many different types of bed frames.



| Article    | Capacity | Tube length | Anti-reflux valve | Article number    | PZN      |
|------------|----------|-------------|-------------------|-------------------|----------|
| UB 2000USo | 2000 ml  | 90 cm       | yes               | 04/07/2000/90/GHC | 08845501 |
| UB 1500USo | 1500 ml  | 90 cm       | yes               | 04/07/1500/90/01  | 16839791 |
| SB 1500USo | 1500 ml  | 90 cm       | no                | 04/07/1500/90/00  | 16839785 |

CARE FLOW UB 1000 Closed urine drainage systems  
Leg bag systems for day & night  
Capacity: 1000 ml

Product specifications:

- ✓ Closed urine drainage system
- ✓ Bacteria- and fluid-tight venting filters provide for permanent pressure compensation and optimal drainage in the system.
- ✓ A doghouse on the bag ensures safe and hygienic positioning of the drain tap.
- ✓ T-drain tap.
- ✓ Needle free sampling port.
- ✓ Clear gradvation marks.



| Capacity | Tube length | Article number    | PZN      |
|----------|-------------|-------------------|----------|
| 1000 ml  | 30 cm       | 04/03/1000/30/GHC | 06718885 |
| 1000 ml  | 60 cm       | 04/03/1000/60/GHC | 06718891 |
| 1000 ml  | 90 cm       | 04/03/1000/90/GHC | 06718916 |

48

UB 2000USo: **Packaging Unit:** 150 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces  
UB 1500USo: **Packaging Unit:** 200 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces  
SB 1500USo: **Packaging Unit:** 200 pieces, **Packaging:** non-sterile **Sub-packaging:** each 10 pieces

Holder for UB 2000USo, UB 1500USo and SB 1500USo

- ✓ Plastic hanger for bedside bags with twin holder. Equipped with two pins taking up single-use bedside bags hanger holes.



| Article number    | PZN      |
|-------------------|----------|
| 04/99/BHL/PVC/GHC | 06718951 |

**Packaging Unit:** 100 pieces **Sub-packaging:** pieces each  
**Packaging:** non-sterile

Pictures may differ from the original

49



Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



Drip chamber

An anti-reflux valve in the drip chamber provides for reliable prevention of urine reflux and highly efficient protection against ascending migration of bacteria.



Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leak-ages caused by improper use of the piston.

**Packaging Unit:** 30 pieces  
**Packaging:** sterile

Pictures may differ from the original

Fixation for leg bags on page 51-56.

CARE FLOW UB 750 + 750US sterile / non-sterile  
1-way chamber leg bag systems for mobile users  
Capacity: 750 ml

Product specifications:

- ✓ The leg bag with 10 cm and 30 cm tube can conveniently be worn on the thigh.
- ✓ The leg bag with 50 cm tube comes with an additional stepped connector in the package. So the tube can be cut to the individual length and the leg bag can conveniently be worn on the thigh or lower leg.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.

| Capacity | Tube length | Article number   | PZN      | Sterility |
|----------|-------------|------------------|----------|-----------|
| 750 ml   | 10 cm       | 04/01/750/10/GHC | 06718862 | sterile   |
| 750 ml   | 30 cm       | 04/01/750/30/GHC | 09772046 | sterile   |
| 750 ml   | 50 cm       | 04/01/750/50/GHC | 06718879 | sterile   |

Packaging Unit: 10 cm: 75 pieces 30 cm: 60 pieces, 50 cm: 50 pieces

| Capacity | Tube length | Article number     | PZN      | Sterility   |
|----------|-------------|--------------------|----------|-------------|
| 750 ml   | 10 cm       | 04/01/750US/10/GHC | 01972716 | non-sterile |
| 750 ml   | 30 cm       | 04/01/750US/30/GHC | 12341427 | non-sterile |
| 750 ml   | 50 cm       | 04/01/750US/50/GHC | 00226543 | non-sterile |

Packaging Unit: 10 cm: 120 pieces 30 cm: 100 pieces, 50 cm: 100 pieces, Sub-Packaging: each 10 pieces



CARE FLOW UB 750F + UB 750UF sterile / non-sterile  
1-way chamber leg bag systems with flip valve for mobile users  
Capacity: 750 ml

Product specifications:

- ✓ Especially users with limited hand dexterity or reduced power can safely and hygienically empty this leg bag, due to our newly developed flip valve.
- ✓ The leg bag with 50 cm tube comes with an additional stepped connector in the package. So the tube can be cut to the individual length and the leg bag can conveniently be worn on the thigh or lower leg.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.

| Capacity | Tube length | Article number    | PZN      | Sterility |
|----------|-------------|-------------------|----------|-----------|
| 750 ml   | 50 cm       | 04/01/750F/50/GHC | 12371500 | sterile   |

Packaging Unit: 50 cm: 50 pieces

| Capacity | Tube length | Article number     | PZN      | Sterility   |
|----------|-------------|--------------------|----------|-------------|
| 750 ml   | 50 cm       | 04/01/750UF/50/GHC | 12371517 | non-sterile |

Packaging Unit: 50 cm: 100 pieces

Packaging: each 10 pieces



50



Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



Doghouse

A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



Drain tap

The special design of the piston in the drain tap ensures safe and easy operation and prevents leakages caused by improper use of the piston.

Packaging: sterile / non-sterile

Fixation for leg bags on page 51-56.

Pictures may differ from the original

51



Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



Doghouse

A doghouse on the bag ensures safe and hygienic positioning of the drain tap.



Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



Drainage tap

The ergonomically formed flat flip valve of the drainage tube can be used easily and intuitively. To drain the leg bag the lever of the valve is flipped downwards by 180°.

Packaging: sterile / non-sterile

Fixation for leg bags on page 51-56.

Pictures may differ from the original

CARE FLOW UB 500 + UB 500US sterile / non-sterile  
3-way chamber leg bag system for mobile users  
Capacity: 500 ml

Product specifications:

- ✓ The leg bag with 10 cm and 30 cm tube can conveniently be worn on the thigh.
- ✓ The leg bag with 50 cm tube comes with an additional stepped connector in the package. So the tube can be cut to the individual length and the leg bag can conveniently be worn on the thigh or lower leg.
- ✓ An anti-reflux valve provides for reliable prevention of urine reflux and efficient protection against ascending migration of bacteria.



| Capacity | Tube length | Article number   | PZN      | Sterility |
|----------|-------------|------------------|----------|-----------|
| 500 ml   | 10 cm       | 04/01/500/10/GHC | 09084990 | steril    |
| 500 ml   | 30 cm       | 04/01/500/30/GHC | 12341203 | steril    |
| 500 ml   | 50 cm       | 04/01/500/50/GHC | 09085009 | steril    |

Packaging unit: 10 cm: 75 pieces, 30 cm: 60 pieces, 50 cm: 50 pieces

| Capacity | Tube length | Article number     | PZN      | Sterility   |
|----------|-------------|--------------------|----------|-------------|
| 500 ml   | 10 cm       | 04/01/500US/10/GHC | 09085015 | non-sterile |
| 500 ml   | 30 cm       | 04/01/500US/30/GHC | 12341410 | non-sterile |
| 500 ml   | 50 cm       | 04/01/500US/50/GHC | 09085021 | non-sterile |

Packaging unit: 10 cm: 120 pieces, 30 cm: 100 pieces, 50 cm: 100 pieces, Sub-Packaging: each 10 pieces

52




Stepped connector

The special design of the connector steps provides for a reliable connection with the ring profile in the funnel end of the catheter. Needle-free urine sampling avoids needle injuries.



3-chambers

Due to the 3-chambers will the urin spread more regular through the system. In that way is it possible to wear the filled bag even more discrete beneath the clothing.



Fleece coating

The backside of the leg bag is coated with sweat-absorbing fleece. This effectively prevents irritation of the skin and considerably improves user comfort.



Drainage tap

The ergonomically formed flat flip valve of the drainage tube can be used easily and intuitively. To drain the leg bag the lever of the valve is flipped downwards by 180°.

Packaging: sterile / non-sterile

Fixation for leg bags on page 51-56.

Pictures may differ from the original

CARE FIX  
Universal Velcro strap  
Length can be individually adjusted

Product specifications:

- ✓ Universal leg bag holding band with two attached Velcro straps to go through the upper eyelets of the leg bag and a special sewn on fixation for the tube of the leg bag system.
- ✓ The elastic material ensures that disturbances of the blood circulation in the thigh and lower leg are prevented as far as possible.
- ✓ The Velcro strap is 70 cm long and can be cut to the individual length. So it can be used on the thigh and lower leg.
- ✓ This universal Velcro strap can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



53

- 60°
- Do not put in the dryer
- Do not iron



| Lenght | Article number    | PZN      |
|--------|-------------------|----------|
| 70 cm  | 05/01/BKL/Uni/GHC | 06718968 |

Packaging Unit: 30 pieces  
Packaging: sterile

Pictures may differ from the original

CARE FIX  
Velcro strap premium set  
Length can be individually adjusted

Product specifications:

- ✓ Sets consists of:  
an extra wide holding band with two sewn-on  
Velcro straps to fixate the upper part of the leg bag.
- ✓ Both sets are equipped sewn-on loop to fixate the  
leg bag tube.
- ✓ The second Velcro Strap ensures a safe fixation of  
the lower part of the leg bag.
- ✓ Silicone strips on the inner side of the holding band  
ensures a good grip on the leg.
- ✓ The elastic material prevents a disruption of the blood  
circulation on the thigh or lower leg.
- ✓ A 70 cm long and individually adjustable Velcro strap can  
be worn on both the thigh as well as the lower leg.
- ✓ The attached loop for an easier opening of the holding  
band.

| Lenght | Article number    | PZN      |
|--------|-------------------|----------|
| 70 cm  | 05/01/BKL/PRE/SET | 09241099 |

**Packaging Unit:** 10 pieces  
**Packaging:** non-sterile



CARE FIX  
Leg pants with integrated leg bag pocket

Product specifications:

- ✓ Soft and ribbed lag pants with integrated  
leg bag pocket.
- ✓ Comfortable and safe fixation on the High  
as well as the lower leg.
- ✓ Due to the special design of the drain tap on  
the leg bag, emptying the bag is very simple.
- ✓ The extra-wide lateral Velcro fastener makes  
it very easy for the user to put on the leg pan  
(Art.-Nr. 05/02/10).
- ✓ The CARE FIX leg pants can be used for the  
following urine leg bag systems:  
UB1000, UB750/UB750US, UB750F/UB750UF,  
UB500/UB500US.

| Size | Article number      | PZN      |
|------|---------------------|----------|
| XXL  | 05/02/10/OS/XXL/GHC | 08880957 |
| XL   | 05/02/10/OS/XL/GHC  | 08871272 |
| L    | 05/02/10/OS/L/GHC   | 08876513 |
| M    | 05/02/10/OS/M/GHC   | 08845748 |
| S    | 05/02/10/US/S/GHC   | 07126649 |

**Measurement:** XXL = 60 - 85 cm, XL = 50 - 60 cm,  
L = 40 - 60 cm, M = 35 - 45 cm, S = 25 - 40 cm

**Packaging Unit:** 10 pieces  
**Packaging:** non-sterile

- 60°
- Do not put in the dryer
- Do not iron

Measure in the middle of the upper or lower leg.

| Size | Article number     | PZN      |
|------|--------------------|----------|
| XL   | 05/02/11/OS/XL/GHC | 12418555 |
| L    | 05/02/11/OS/L/GHC  | 07126017 |
| M    | 05/02/11/OS/M/GHC  | 07126023 |
| S    | 05/02/11/US/S/GHC  | 07126313 |

**Measurement:** XL = 65 - 95 cm, L = 55 - 80 cm,  
M = 45 - 65 cm, S = 30 - 45 cm

**Packaging Unit:** 10 pieces  
**Packaging:** non-sterile



54

- ✓ This Velcro strap premium set can be used for the  
following urine leg bag systems: UB1000, UB750/  
UB750US, UB750F/UB750UF, UB500/UB500US.
- ✓ The premium set with buttons (Art.-Nr. 05/01/BKL/PRE/03)  
comes with additional four buttons to fixate the upper  
part of the leg bag.

- 60°
- Do not put in the dryer
- Do not iron

| Lenght | Article number   | PZN      |
|--------|------------------|----------|
| 70 cm  | 05/01/BKL/PRE/03 | 12357931 |

**Packaging Unit:** 10 pieces  
**Packaging:** non-sterile



55

Pictures may differ from the original

Pictures may differ from the original

CARE FIX  
Thigh cuffs sizes: XL, L, M

Product specifications:

- ✓ The fabric leg bag holder made from a washable cotton blend has an integrated leg bag pocket.
- ✓ A silicon strip on the inner side ensures a good grip.
- ✓ The fabric leg bag holder can be worn on both the right as well as the left thigh.
- ✓ The bags are positioned into the integrated pocket with a Velcro fastener; including an opening for the drain tap.
- ✓ The CARE FIX cuffs can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



CARE FIX  
Lower leg cuffs sizes: L, M, S

Product specifications:

- ✓ The fabric lower leg bag holder made from a washable cotton blend has an integrated leg bag pocket.
- ✓ A silicon strip on the inner side ensures a good grip.
- ✓ The fabric lower leg bag holder can be worn on both the right as well as the left lower leg.
- ✓ Leg bags are positioned into the integrated pocket with a Velcro fastener; including an opening for the drain tap.
- ✓ The CARE FIX cuffs can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



56

- 60°
- Do not put in the dryer
- Do not iron

It is measured in the middle of the thigh.

Leg bag holder for the thigh

| Size | Article number  | PZN      |
|------|-----------------|----------|
| XL   | 05/02/12/XL/GHC | 09080408 |
| L    | 05/02/12/L/GHC  | 09080383 |
| M    | 05/02/12/M/GHC  | 09080377 |

Measurement: XL = 61 cm, L = 58 cm, M = 54 cm

Packaging Unit: 10 pieces

Packaging: non-sterile



Pictures may differ from the original

57

- 60°
- Do not put in the dryer
- Do not iron

It is measured in the middle of the lower leg.

Lower Leg Cuffs

| Size | Article number | PZN      |
|------|----------------|----------|
| L    | 05/02/13/L/GHC | 09080437 |
| M    | 05/02/13/M/GHC | 09080420 |
| S    | 05/02/13/S/GHC | 09080414 |

Measurement: L = 37 cm, M = 35 cm, S = 32 cm

Packaging Unit: 10 pieces

Packaging: non-sterile



Pictures may differ from the original

CARE FIX  
Premium one leg pants sizes: XL, L, M

Product specifications:

- ✓ The premium one leg pants are made from a washable cotton blend has an integrated leg bag pocket.
- ✓ Individually adjustable belt with silicone strips for a safe fixation.
- ✓ The premium one leg pant can be worn on both the right as well as the left thigh.
- ✓ Leg bags are positioned into the integrated pocket with a Velcro fastener; including an opening for the drain tap.
- ✓ Easy emptying of the leg bags due to the comfortable positioning of the drain tap.
- ✓ The leg pants can be used for the following urine leg bag systems: UB1000, UB750/UB750US, UB750F/UB750UF, UB500/UB500US.



CARE FLOW und CARE VENT catheter valve  
For an independent and controlled emptying of the urinary bladder

Product specifications:

- ✓ The optimal solution for mobile, independant users wearing a transurethral or suprapubic catheter.
- ✓ The catheter valve continuously trains the bladder muscles so that the bladder function remains completely intact. The emptying of the urinary bladder follows the user's personal rhythm. This sustainably enhances the quality of life.
- ✓ During the night the valve can be connected to a urine bag system by means of a firmly attached silicone tube.
- ✓ By flipping the lever of the flip valve by 180° the CARE FLOW catheter valve is completely opened as well as fixed in a safe position. This prevents an unintentional closing of the valve during the night.



Arrows indicate open / closed position



Closed catheter valve



Open catheter valve, flipped by 180°

CARE FLOW catheter valve

| Weight | Article number  | PZN      |
|--------|-----------------|----------|
| 3g     | 03/02/10/00/GHC | 12378382 |

Packaging Unit: 5 pieces  
Packaging: sterile

- ✓ The special construction prevents the after-dripping of urine, due to the which ensures a hygienical operation and utmost discretion.
- ✓ A colour coding of the pivot and slide mechanism makes it very easy to use single handed.

CARE VENT catheter valve

| Weight | Article number  | PZN      |
|--------|-----------------|----------|
| 6g     | 03/02/10/01/GHC | 01914676 |

Packaging Unit: 10 pieces  
Packaging: non-sterile



- 60°
- Do not put in the dryer
- Do not iron

Measure in the middle of the upper or lower leg.

| Size | Article number  | PZN      |
|------|-----------------|----------|
| XL   | 05/03/10/XL/GHC | 09080472 |
| L    | 05/03/10/L/GHC  | 09080466 |
| M    | 05/03/10/M/GHC  | 09080443 |

Measurement: XL = 98 - 120 cm, L = 90 - 100 cm, M = 85 - 95 cm  
Packaging Unit: 10 pieces  
Packaging: non-sterile



Pictures may differ from the original

Pictures may differ from the original

CARE DRESS  
Suprapubic dressing set

Product specifications:

- ✓ Suprapubic dressing set for optimal care of the suprapubic exit site.
- ✓ The included smaller plaster serves to temporarily fix the suprapubic catheter to the abdominal wall, thus preventing the catheter from shifting.



CARE DRESS  
Abdominal PEG dressing set

Product specifications:

- ✓ Abdominal PEG dressing set for optimal care of the PEG exit site.



Included in this set:

- 1 x elastic plaster 10 x 10 cm
- 1 x elastic plaster 5 x 3 cm
- 3 x non-woven swabs 5 x 5 cm, 6-ply
- 2 x non-woven split swabs 5 x 5 cm, 6-ply
- Packed in a peel-apart pouch -

| Article number | PZN      |
|----------------|----------|
| 08/02/VS/GHC   | 06718974 |

Packaging Unit: 25 sets  
Packaging: sterile

Pictures may differ from the original

Included in this set:

- 1 x elastic plaster 12 x 15 cm
- 3 x non-woven swabs 7,5 x 7,5 cm, 8 ply
- 2 x non-woven split swabs 7,5 x 7,5 cm, 4 ply
- Packed in a peel-apart pouch -

| Article number | PZN      |
|----------------|----------|
| 08/02/PEG/GHC  | 05461059 |

Packaging Unit: 15 sets  
Packaging: sterile

Pictures may differ from the original



Visit our website [www.ghc-medical.de](http://www.ghc-medical.de).

There you will find detailed information about our company as well as the entire product portfolio.

Our sales team is looking forward to your call:

Mon.-Thurs.: 08:00 – 17:00  
Fri.: 08:00 – 14:00

**Contact:**

GHC German Health Care GmbH  
Birkstück 15  
24999 Wees/Flensburg

Tel.: +49 (0) 46 31 / 57 200-0  
Fax: +49 (0) 46 31 / 57 200-19

[info@ghc-medical.de](mailto:info@ghc-medical.de)  
[www.ghc-medical.de](http://www.ghc-medical.de)